RANGE OF MOTION ADJUSTMENTS

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ENGLISH

The Stride OA comes with standard 5 degree extension stops. Additional flexion and extension stops may be installed if they are prescribed by your physician or needed for optimal performance or fit of the brace. These should be inserted by the technician who fits your brace.

Optional 0°, 5°, 10°, 15°, 20°, and 25° extension stops and optional 45°, 60°, 75° and 90° flexion stops are included with each brace.

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CALF SUSPENSION WRAP APPLICATION INSTRUCTIONS

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ENGLISH

A suspension wrap may have been provided with your brace if necessary. This wrap is designed to aid suspension and brace position by applying even compression to soft tissue.

STEP 1.

Attach self-adhesive Velcro strip to calf cuff and dot to lower upright.

STEP 2.

Attach wrap end to Velcro on cuff and lower upright (neoprene side faces skin) and pull snugly around back of leg.

STEP 3.

Stretch the wrap around the leg and attach it to itself.

NOTE: Wrap should be positioned so that top edge of wrap is 1/2" inch above top edge of calf strap #1.

CARE AND MAINTENANCE

Wrap should be removed and hand-washed as needed using cold water and mild detergent. Rinse thoroughly and air-dry (do not machine-dry).

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ENGLISH

INDICATIONS:

- Mild to severe unicompartmental knee osteoarthritis
- Medial or lateral compartment

WARNINGS:

If you experience pain, swelling, sensation changes, or unusual reactions while using this product, consult a physician.

FOR SINGLE PATIENT USE ONLY.

NOTICE: While every effort has been made using state-of-the-art techniques to obtain maximum compatibility of function, strength, durability and comfort, this device is only one element in the overall treatment program administered by a medical professional. There is no guarantee that injury will be prevented through use of this device.

CAUTION: Applicable federal, state, and/or local regulations may restrict this device to sale by or on the order of a practitioner licensed by law of the State in which he/she practices to use or order the use of this device.

WARRANTY: Corflex warranties all its products from the original date of purchase against defects in materials and workmanship. Normal wear and tear during use of a product is not considered a defect. Contact Corflex at **800.426.7353** for specific product warranty information.

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STEP 2 / PASO 2 / ÉTAPE 2

STEP 3 / PASO 3 / ÉTAPE 3

STRIDE OA APPLICATION INSTRUCTIONS

(Instrucciones de Aplicación de Stride OA / Instructions de Mise en Place de Stride OA)

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ENGLISH

STEP 1. LOOSEN STRAPS 1 - 4

- To simplify application, attach each strap back onto itself to minimize tangling.

STEP 2. POSITION BRACE

- While sitting at edge of chair, bend knee slightly and apply brace to lateral side (outside) of leg.
- Align center of hinge with top of kneecap and position slightly behind midline of leg.

Important note for initial fitting: When positioning the **Stride OA off-the-shelf brace**, it should first be adjusted to neutral. This is achieved by setting the correction angle (see step 5) so the cuffs and hinge rest on leg without gapping or pressure. Medial or lateral compartment correction will be adjusted in step 5.

STEP 3. FASTEN STRAPS

Attach straps in sequence indicated on strap ends. Connect buckles then tighten straps. Hold brace at hinge to prevent rotation while tightening straps.

- Fasten strap 1 snugly just above calf muscle. Strap 1 is important in maintaining brace position.
- Fasten strap 2 around middle of calf.
- Fasten strap 3 around lower thigh.
- Fasten strap 4 around upper thigh.

All straps should be comfortably tightened to properly secure brace and provide relief.

NOTE: Elastic strap extension may be used on strap 3 for medial compartment and strap 4 for lateral compartment. Remove strap 3 or 4. Place elastic extension through frame slot and attach to strap as shown.

STEP 4. CHECK FIT

- Once straps are secure, center of hinge should be aligned with middle of kneecap and positioned slightly behind midline of leg.
- To ensure proper fit for activity, walk a short distance. It is normal for brace to drop or settle into position. Make strap and position adjustments as necessary.

STEP 5. ADJUST BRACE CORRECTION ANGLE (Q-HINGE)

Warning: Initial adjustment of the Stride OA brace must be done by a trained bracing professional in accordance with a Doctor's prescription indicating the compartment to unload. Subsequent adjustments must be made only by a trained bracing professional or by the patient if the patient has been given detailed instruction from a professional and has fully read and understands these application instructions.

For initial fitting of Stride OA off-the-shelf: With the brace in neutral (from step 2 above), bend the knee to 90°. With the provided allen key, turn the Q-hinge screw 1/8 turn clockwise (**important**) to increase lateral compartment relief or 1/8 turn counter-clockwise to increase medial compartment relief.

This is a good starting point for initial use.

STEP 5.1 FOR INITIAL FITTING AND ALL ADJUSTMENTS

- Stand and walk a short distance to ensure comfort and proper correction after each adjustment. If there is any discomfort in the knee joint or leg, turn the angle adjustment back slightly.
- If more relief is required, turn Q-hinge screw no more than 1/8 turn. Walk after the adjustment and check as above.
- Over-adjusting the Q-hinge may cause discomfort and skin pressure issues. Adjust conservatively.

“WEAR TIME”

Gradually increase brace wear time to allow your body to adjust to the brace. Start with 30 minutes on day one and gradually increase. Pain caused by osteoarthritis may take time to subside. It is important that you routinely wear the brace during weight bearing activities (e.g. when you're on your feet) as it will reduce joint loading and allow you to begin feeling relief. It may take a number of weeks for relief to become significant.

CARE AND MAINTENANCE

Salt water, sand, dirt and other elements can cause damage. Thoroughly rinse brace with fresh water after exposure to these conditions. Towel or air-dry. Liners should be removed and hand-washed as needed using cold water and mild detergent. Rinse thoroughly and air-dry (do not machine-dry.) Lubricate hinges with dry silicone spray as needed.

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© 2014 Corflex, Inc.