



CORFLEX®

Designed for life.

CORFLEXGLOBAL.COM | 800 426 7353

CoolTex™ Ag Hybrid Anterior Closure Knee with ROM Hinge (13" & 16") Application Instructions

Fitting and Application:

1. Open thigh closure and unfasten hook and loop straps. Determine proper positioning of Cor-Trak buttress depending on patient need (inferior, superior or lateral).



2. If necessary remove hinge from brace and adjust range of motion stops (comes preset at 10° of extension). Both hinges must have the same settings for proper function.



3. Slide foot through the sleeve and locate bottom of patella opening with bottom of the knee cap. Make sure center axis of hinge is aligned with the center of the knee joint.



4. Secure closure across top of thigh to desired comfort level.



5. Starting with strap nearest the foot, secure strap through d-ring. Tighten to desired comfort level. Attach remaining straps in the same manner.



Care Instructions

The CoolTex Ag Hybrid Anterior Closure Knee with ROM Hinge may be hand washed using warm water and detergent for delicate fabrics. Remove ROM Hinge prior to washing. Repeatedly rinse and line dry. Do not use bleach or place in a washer or dryer.

Caution

This product has been designed and tested based on single patient use. Since it is a single use device it should never be used by multiple patients. If any additional pain or symptoms occur while using this device, seek medical attention.

Indications

Indicated for medial/lateral instability, muscle/ligament weakness or injury, ACL/PCL moderate instability, strains or sprains, meniscal injuries and degenerative joint disease.

Contraindications

Not for use on patients with circulatory or sensation issues. If swelling, pain, skin irritation or an unusual reaction occurs, discontinue use and consult your medical professional.

Warranty

Products manufactured by Corflex are warranted to the original purchaser against defects in workmanship and material for the life of the product. The foregoing warranty shall not apply to normal wear and tear; or if the product has been subjected to misuse or misapplication the warranty will be voided at the discretion of Corflex Global.

To Reorder

Product Number

	XS	S	M	L	XL	2XL	3XL	4XL
13"	85-8322	85-8323	85-8324	85-8325	85-8326	85-8327	85-8328	85-8329
16"	85-8342	85-8343	85-8344	85-8345	85-8346	85-8347	85-8348	85-8349

Sizing Guide

6" Above Patella	13"-16"	14"-18"	15"-20"	16"-23"	17"-25"	18"-26½"	19"-28½"	21"-30½"
Mid Patella	11"-12"	12"-13"	13"-14"	14"-16"	16"-18"	18"-20"	20"-22"	22"-24"
6" Below Patella	11"-12"	12"-13"	13"-14"	14"-15"	15"-16"	16"-17"	17"-18"	18"-19"

Accessories

Cor-Trak Buttress 10/Pk	88-9000
Cor-Trak Buttress	88-9001



669 EAST INDUSTRIAL PARK DRIVE
 MANCHESTER, NH 03109
 CORFLEXGLOBAL.COM | 800 426 7353