

Marathon Speed Lacer™ Ankle Stabilizer

Superior compression supports performance and recovery

Marathon Speed Lacer ankle brace with speed lacer closure allows fast and easy application while the stabilizing ankle straps offer support and protection. The low-profile design fits any type of shoe and the CoolMesh Ag™ enhances breathability and long-term wear.



TexSmart™

Featuring one of our technology-based materials, CoolMesh Ag™, the Marathon Speed Lacer Ankle Stabilizer delivers innovative, high performance solutions for a wide range of uses.



CORFLEX®


Designed for life.



Marathon Speed Lacer™ Ankle Stabilizer

Made with high-quality ballistic nylon material, a thin yet durable nylon weave that maintains its shape while providing maximum support.

- Speed Lacer closure allows quick application and adjustment of the brace.
- Hook and loop closure reduces the chance of accidental loosening for a secure fit.
- Features dual nylon figure-8 stabilizing straps that replicate ankle taping and lock the calcaneus to control abnormal inversion/eversion.
- CoolMesh Ag's ventilated, antimicrobial fabric tongue and Achilles Comfort panel provide comfort and protection; wide elastic Tib/Fib strap for added proximal compression.
- Indicated for chronic ankle instability, sprains or strains and protection against inversion or eversion ankle injuries.
- Constructed of ballistic nylon, CoolMesh Ag, HDPE removable stays, polyester webbing, latex-free elastic.



CoolMesh Ag™
Freshness is in the Silver

Breathable CoolMesh encourages airflow and wicks moisture away from body up to 6 times faster than other fabrics.

CoolMesh Ag's specialty silver yarn permanently inhibits the growth of odor-causing bacteria, mold, mildew and algae on fabric.

The addition of the antimicrobial silver helps boost the fabrics durability and lifespan.

Customized stabilization with removable padded and contoured stays



Marathon Speed Lacer Ankle Stabilizer

Product Number

2XS	XS	S	M	LG	XLG	2XLG	3XLG
75-8130	75-8131	75-8132	75-8133	75-8134	75-8135	75-8136	75-8137

Sizing Guide Circumference

9"-10"	10"-11"	11"-12"	12"-13"	13"-14"	14"-15"	15"-16"	16"+
--------	---------	---------	---------	---------	---------	---------	------

Measurement is circumference from base of heel around ankle.
Fits either left or right foot.

Sizing Guide - Shoe Size

Men	6-7	7½-9	9½-11	11½-13	13½-15	15+
Narrow	2XS	XS	S	M	L	XL
Regular	XS	S	M	L	XL	2XL
Wide	S	M	L	XL	2XL	3XL
Women	5-6	6½-8	8½-10	10½-12	12½-14	14+
Narrow	2XS	XS	S	M	L	XL
Regular	XS	S	M	M	L	XL
Wide	S	M	L	LL	XL	2XL

Suggested Billing Code L-1902 Ankle Foot Orthosis



All products manufactured by Corflex are Latex Free.

CORFLEX | 669 EAST INDUSTRIAL PARK DRIVE | MANCHESTER, NH 03109
CORFLEX.COM | TOLL FREE: 800 426 7353 | FAX: 603 623 4111

© 2020 Corflex, Inc. All Rights Reserved. Printed in USA. 1360/6-20



CORFLEX®
Designed for life.