

# Rival Pro and Rival Ace

# Rival Pro and Rival Ace Ligament Knee Braces

Ready to take on your toughest competitor



Cotton/Lycra Undersleeve



Cotton/Mesh Oversleeve

## Rival Pro Series

Product Number	S	M	L	XL	2XL
ATC Top Right	81-0401	81-0402	81-0403	81-0404	81-0405
ATC Top Left	81-0411	81-0412	81-0413	81-0414	81-0415
Rigid Top Right	81-1421	81-1422	81-1423	81-1424	81-1425
Rigid Top Left	81-1431	81-1432	81-1433	81-1434	81-1435

Accessories	S	M	L	XL	2XL
Cl Strap 20"L	81-0630				
Strap Kit	81-0451	81-0452	81-0453	81-0454	81-0455
Liner Kit ATC Right	81-0461	81-0462	81-0463	81-0464	81-0465
Liner Kit ATC Left	81-0471	81-0472	81-0473	81-0474	81-0475
Liner Kit Rigid Top Right	81-1481	81-1482	81-1483	81-1484	81-1485
Liner Kit Rigid Top Left	81-1491	81-1492	81-1493	81-1494	81-1495

## Rival Ace Series

Product Number	S	M	L	XL	2XL
ATC Top Right	81-0501	81-0502	81-0503	81-0504	81-0505
ATC Top Left	81-0511	81-0512	81-0513	81-0514	81-0515
Rigid Top Right	81-1521	81-1522	81-1523	81-1524	81-1525
Rigid Top Left	81-1531	81-1532	81-1533	81-1534	81-1535

Accessories	S	M	L	XL	2XL
Cl Strap 20"L	81-0630				
Strap Kit	81-0551	81-0552	81-0553	81-0554	81-0555
Liner Kit ATC Right	81-0561	81-0562	81-0563	81-0564	81-0565
Liner Kit ATC Left	81-0571	81-0572	81-0573	81-0574	81-0575
Liner Kit Rigid Top Right	81-1581	81-1582	81-1583	81-1584	81-1585
Liner Kit Rigid Top Left	81-1591	81-1592	81-1593	81-1594	81-1595

## Sizing Guide

6" Above Patella	14½"-17"	17"-19½"	19½"-22½"	22½"-25½"	25½"-30"
Mid Patella	11½"-13"	13"-15"	15"-17"	17"-19"	19"-22"
6" Below Patella	11½"-13"	13"-14½"	14½"-16½"	16½"-18½"	18½"-22"

Measurements may fall into different size ranges. If patient is between two sizes, larger model is typically most appropriate.

## Suggested Billing Code

L-1845 (Custom Fit) or L-1852 (OTS)  
Knee Orthosis Double Upright, Prefabricated

## Sleeves

Undersleeve	OverSleeve	
XS	80-0010	
S	80-0011	80-0031
M	80-0012	80-0032
L	80-0013	80-0033
XL	80-0014	80-0034
2XL	80-0015	80-0035
3XL	80-0016	80-0036
4XL	80-0017	80-0037

## Sizing Guide

XS	13"-20"	
S	14"-22"	Up to 15"
M	15"-24"	15"-17¾"
L	16"-26"	17¾"-19½"
XL	19"-28"	19½"-21½"
2XL	22"-29"	21½"-24½"
3XL		24½"-28"
4XL		28"-34"

Circumference measurement taken 8" above mid-patella



All products manufactured by Corflex are Latex Free.



Rival Ligament Braces feature four styles to best match your patient's competitive needs. Whether mild to severe ACL, PCL or combined instabilities, Rival Ligament Braces deliver solid support and stability, while providing superior protection for contact and non-contact sports and activities.



**CORFLEX**  
Designed for life.



# Rival Pro and Rival Ace

Ligament braces are indicated for ACL, PCL, MCL, LCL or combined ligament instabilities. Used post-surgically to protect repairs or as conservative treatment of ligament damage or deficiency to help provide stability, support and protection during rotational and anteroposterior activities.



Rival Pro Ligament Brace with Rigid Top Cuff



Rival Pro Ligament Brace with Active Thigh Cuff

## Rival Pro

The Rival Pro Ligament Brace is designed with a posterior calf cuff for improved tibial crest comfort and maximum protection. The posterior cuff along with anterior tibial straps control anterior tibial translation while maintaining a high level of comfort by minimizing tibial pressure.



Rival Ace Ligament Brace with Rigid Top Cuff



Rival Ace Ligament Brace with Active Thigh Cuff

## Rival Ace

The Rival Ace Ligament Brace features an anterior cuff to comfortably fit the tibial crest while providing superior rotational control. The anterior tibial cuff with posterior straps allows muscles to move freely minimizing brace migration and pistoning.



### Active Thigh Cuff

The patented Active Thigh Cuff™ (ATC) dynamically adapts to the shape and contour of the thigh throughout the full range of motion avoiding pistoning and migration seen in traditional designs. Flexible straps easily adjust allowing for variable sizing and enhance brace suspension, protection and comfort.



### Rigid Top Cuff

The Rigid Top Cuff helps deliver maximum stability to torn or injured ligaments and offers additional stability and exceptional protection during contact sports in an easy to apply, low profile brace. Both models feature a shorter medial thigh upright for improved comfort.

## Delivering winning results

Superior suspension and comfort without compromising protection. Enables a customized fit that moves with the muscles through the full range of motion and avoids the pistoning and migration seen in traditional designs. Along with the Rival Ace's anterior tibial cuff, you will achieve flexibility of movement while maintaining control and stability.

### Enhanced Strapping System

Strapping system incorporates swivel attachments that contour to the leg to help enhance suspension and comfort. Molded strap tips allow for easy removal and application.



### Active Thigh Cuff

Adapts to the shape and contour of the thigh throughout the full range of motion.

### Superior Brace Suspension

Compression molded faux-suede padded liners and specialty strapping feature silicone for excellent brace suspension.



### Adjustable Knee Extension

Metal knee hinges feature movable stops that control knee extension from 0° to 25° in 5° increments and flexion from 45° to 90° in 5° increments.

### Built-in Durability

Strong yet low-profile, aerospace-grade aluminum uprights are 30% thicker than other aluminum braces.

### High Compression-Set Resistance

Lower liner cushions conform to the tibial crest and bounce back from compression even under extreme conditions.



### Rigid Anterior Tibial Cuff

Provides maximum control of tibial rotation and anterior drawer.



Rival Ace with Active Thigh Cuff