



CORFLEX®

Designed for life.

800 426 7353

Cryotherm™ Hip Wrap

Our latest addition to the Cryotherm line of economical cold/hot therapy which includes knee, shoulder, elbow, wrist, ankle, back and neck wraps.

CORFLEX.COM



Cryotherm™ Hip Wrap

Anatomically contoured wrap provides mild compression and concentrated cold/hot therapy to the hip for post-surgical and non-surgical indications. Soft breathable foam laminate wrap incorporates two interior cloth pockets that hold the Cryotherm reusable gel packs securely in place and offer a protective barrier between the skin and gel pack.

Wide waist strap makes application easy and adjustable hip closures allow for a customized fit. Wrap provides a broad coverage area (11" x 22") and can be rotated for use with anterior, lateral or posterior approaches or conditions. Universal sizing fits left or right to minimize inventory requirements.

- > Designed to contour to the hip and provide mild compression with the benefits of cold/hot therapy; wrap can be positioned to provide coverage to the anterior, lateral or posterior hip anatomy
- > Indicated post-surgically following Total Hip Replacement, Resurfacing, Revision or Arthroscopy; non-surgically following injections to the hip joint, arthritis, bursitis, tendonitis, muscle or tendon strains, labral tears, femoroacetabular impingement, piriformis syndrome, sciatica and hip dislocation/subluxation
- > Manufactured from a soft breathable foam laminate with contact closures and pockets to hold two reusable gel packs



Product Number	0 Gels	2 Gels	4 Gels
	72-7060	72-7062	72-7064

Sizing Guide

Sizing Guide
Universal
Fits up to 54" waist circumference and 36" thigh circumference.

Accessories

Gel Pack, 5 1/2" x 11", each	72-7001
Gel Pack, 5 1/2" x 11", 10/bx	72-7000

Visit www.corflex.com/coldtherapy for our full line of Cryotherm and Cryo Pneumatic Rehabilitation Products



PRODUCTS ARE LATEX FREE



PRECAUTIONS NEEDED WHEN USING THIS PRODUCT