

Ultra Stretch Plantar Fasciitis Night Splint



CORFLEX®

Designed for life.

CORFLEX.COM

- Features bi-lateral dorsiflexion assist straps that provide a controlled amount of dorsiflexion and removable wedge for added toe flexion; slip resistant tread
- Indicated for night use of treatment of plantar fasciitis, Achilles tendonitis, plantar flexion contractures, heel spurs and other lower leg over-use injuries
- Constructed of polyethylene shell with foam liner and padded straps with quick-release buckles; textured vinyl tread
- Fits both left or right foot

ADDITIONAL FEATURES

- QUICK RELEASE BUCKLES MAKE APPLICATION AND REMOVAL EASY
- PADDED ADJUSTABLE STRAPS
- BI-LATERAL STRETCH STRAPS FOR CUSTOM DORSIFLEXION
- FOAM TOE WEDGE FOR ADDED PLANTAR FASCIA STRETCH
- BOTTOM TREAD PROMOTES SAFETY
- THIS PRODUCT IS NOT MADE FROM NATURAL RUBBER LATEX



	XS	S	M	L	XL
Product Number					
	37-4330	37-4331	37-4332	37-4333	37-4334
Sizing Guide					
Mens Size	3 - 4 1/2	5 - 6 1/2	7 - 9	9 1/2 - 11	11 1/2 - 15
Womens Size	4 - 5 1/2	6 - 7 1/2	8 - 10	10 1/2 - 12	12 1/2 - 16

Measure By US Shoe Size



Products are Latex Free



CORFLEX®
Designed for life.

CORFLEX | 669 EAST INDUSTRIAL PARK DRIVE | MANCHESTER, NH 03109
CORFLEX.COM | TOLL FREE: 800 426 7353 | FAX: 603 623 4111