



CORFLEX®

Designed for life.

800 426 7353

CORFLEX.COM

Universal Wrist Lacer Orthoses

Reduce your inventory by 75%. Choose right or left, then fit to most patients in just seconds with the easy-to-use cinch closure.



Universal Lacer Wrist



Universal Lacer Wrist w/Abducted Thumb

Universal Lacer Wrist and Wrist w/Abducted Thumb Orthoses are ideal for offices with limited space. Inventory is reduced from 10 items down to two—right and left. Hook and loop dorsal stay assembly can be repositioned to fit a wide range of patients. Lacer closure and stockinette interior allow for easy slip-on application and removal. Malleable palmar stay is contoured and removable. Wrist/Thumb model features two additional malleable stays for greater immobilization.

Universal Lacer Wrist Orthosis

- > Anatomically designed and made from high quality materials for enhanced durability and patient compliance; universal sizing reduces inventory
- > Indicated for Carpal Tunnel Syndrome, tendonitis, sprains or strains, soft tissue injury, wrist contracture and post-cast removal
- > Constructed of breathable, cottony soft, moisture-wicking liner and high performance brushed nylon laminate exterior; malleable dorsal and palmar aluminum stays

Universal Lacer Wrist w/Abducted Thumb Orthosis

- > Adjustable cinch closure features two separate pull straps for a secure fit with minimal effort; malleable thumb stay provides extension and abduction control; contoured palmar stay supports the wrist and is easily modified
- > Indicated for deQuervain's Syndrome, tendonitis, ulnar collateral ligament injury, CMC arthritis, scaphoid injuries, sprains or strains and post-cast removal
- > Constructed of breathable, cottony soft, moisture-wicking liner and high performance brushed nylon laminate exterior; four malleable aluminum stays

Product Number	Right	Left
Wrist Orthosis	73-7500	73-7501
Wrist Orthosis w/ Abducted Thumb	73-8500	73-8501

Sizing Guide

Fits patients with a wrist circumference measuring 5" to 9"

Suggested Billing Codes:

Wrist Orthosis: L-3908 Wrist Hand Orthosis, OTS

Wrist Orthosis w/Abducted Thumb: Wrist Hand Finger Orthosis L-3807 (Custom Fit)
L-3809 (OTS)