

Contender™ Post-Op Elbow Brace

Built to limit motion, not life.



Optional Cradle Arm Sling (included)



Hand Attachment (sold separately)

Indicated for use following elbow reconstruction, ligament and tendon repairs (Tommy John surgery, distal biceps or triceps tendon repair, etc.), collateral ligament reconstructions, elbow hyperextension, tennis elbow release, stable or internal fixed fractures of the distal humerus or proximal to middle radius or ulna or mild to severe soft tissue injuries.

- > **17 oz.** The lightest Post-Op Elbow Brace featuring conformable metal cuffs
- > **Removable Biceps Cuff** Available only on the Contender, helps avoid incision site irritation
- > **QuikZip™ Buckle System** Allows one-handed patient application, making it the easiest post-op elbow brace to apply
- > **ZipTrak™ System** Facilitates length adjustment from 14"-19.5" in just seconds for a precise fit

Product Number	Right	Left
	75-4000	75-4001

Overall length is adjustable from 14" to 19½". Fits biceps circumference up to 22" and forearm circumference up to 18".

Accessories	Universal	Right	Left
Replacement Sling	75-4032		
Replacement Pad Kit		75-4030	75-4031

Hand Attachment	Right	Left
Small	75-4011	75-4021
Medium	75-4012	75-4022
Large	75-4013	75-4023

Hand Attachment Sizing Guide

Small	6"-7"
Medium	7"-8"
Large	8"-9"

Measurement taken from wrist circumference

Suggested Billing Code L-3760 Elbow Orthosis with Joint



CORFLEX®

Designed for life.

800 426 7353 / CORFLEX.COM

669 EAST INDUSTRIAL PARK DRIVE / MANCHESTER, NH 03109

FAX 603 623 4111 / EMAIL SERVICE@CORFLEX.COM



ALL PRODUCTS MANUFACTURED BY CORFLEX ARE LATEX FREE.



CORFLEX®

Designed for life.

Contender™ Post-Op Elbow Brace

Innovative features like the removable biceps cuff and easy, one-handed application make the Contender the clear winner for controlled ROM at the elbow.



Contender™ Post-Op Elbow Brace Built to limit motion, not life.

ZipTrak™ Lock

Allows quick, tool-free adjustment of brace length from 14" to 19.5"

Removable Biceps Cuff

Helps avoid incision site irritation and promotes greater flexion range



QuikZip™ Buckles

Allow for easy application and micro adjustment of strap for a comfortable and secure fit

Trim to Fit Straps

Simply cut to length and apply alligator-style hook to ends

ZipTrak™ System

Adjust overall length as well as individual cuff positions in 1/8" increments to avoid surgical sites

Drop Lock Button and Degree Indicator

11 locking positions to quickly immobilize the brace from -10° to 90°

Ulnar Styloid Relief

Pre-shaped relief helps eliminate hot spots at the ulnar styloid

FlexLok™ & ROM Buttons

Finger touch controls adjust range of motion from -10° to 110° and feature a unique FlexLok band to help prevent patient tampering

Metal Uprights

Custom bendable, aircraft aluminum uprights are notched to assist in contouring

Conforming Cuffs

Malleable aluminum cuffs are easily shaped to the arm to limit brace migration

Compression Molded Padding

Removable fabric covered pads feature comfort channels that encourage airflow and help wick moisture



CORFLEX®

Patent pending.