

Marathon™ Air Walkers

Lightweight Support - Low Profile Design

A circumferential shell with pneumatic liner combines for customized stabilization and edema control.



CORFLEX®

Designed for life.

Marathon™ Air Walkers

Lightweight Support - Low Profile Design

The Marathon Air Walker features the latest innovative design in lightweight comfort with a durable, semi-rigid shell that protects and supports the limb. The air cell in the liner adjusts easily using the built-in inflation bulb and knob. This fully integrated pneumatic air bladder combines with circumferential protection to provide optimum immobilization and edema control.

> Indicated for stable foot fractures, severe sprains, soft tissue injuries, post-op and post-cast removal.

Proportional Scaling

The height and width scale proportionally as boot increases in size for a better fit and improved performance.



Shock Absorbent Insole

A cushioned shock absorbent insole offers comfort in a low-profile design that minimizes leg length discrepancy and foot fatigue.



Low-Profile Rocker Sole

Its low-profile rocker sole is engineered to promote a natural gait while delivering superior cushioning from heel to toe.

Fully integrated Pneumatic Control

An integrated pump and interior bladder allow for a custom fit and increased stabilization throughout recovery.




Durability and Comfort

Features durable, semi-rigid ventilated panels and moisture wicking liner for long-wearing comfort and dryness.



Additional Ankle Strap

Additional ankle strap ensures proper positioning and supplemental support.



D-Ring Technology

D-Rings are attached to the semi-rigid frame making the boot easier to apply and remove.



Removable Toe Guard

A removable and adjustable toe guard provides protection against hazards.



Marathon™ Air Walkers

Suggested Billing Codes
Pneumatic walking boot
L-4360 (Custom Fit)
L-4361 (OTS)



Ankle

Tall



Walker Socks

This thick mid-calf sock is made from 100% synthetic knit and comfortably fits your foot and leg. A non-constricting, seamless design with moisture wicking properties helps prevent irritation. Available separately, 14" H and 18" H.

Suggested Billing Code
Tibial length fracture sock
L-2840



All products manufactured by Corflex are Latex Free.

Marathon Air Walkers - Ankle

Product Number

	XS	S	M	L	XL
Air Walker Ankle	79-2020	79-2021	79-2022	79-2023	79-2024
Replacement Liner Ankle	79-2040	79-2041	79-2042	79-2043	79-2044

UNIV

Ankle Walker Sock (Pair) 79-0274

Sizing Guide *Measured by Shoe Size*

	XS	S	M	L	XL
Men's Size	3-4	4½-7	7½-10	10½-12½	13+
Women's Size	4-5½	6-8	8½-11	11½-13	14+

Marathon Air Walkers - Tall

Product Number

	XS	S	M	L	XL
Air Walker Tall	79-2010	79-2011	79-2012	79-2013	79-2014
Replacement Liner Tall	79-2030	79-2031	79-2032	79-2033	79-2034

UNIV

Tall Walker Sock (Pair) 79-0275

