

CoolTex™ Ag 13" Anterior Closure Knee Wrap with Hinge

Silver-based antimicrobial lining inhibits microbes that can break down fabrics making this the perfect brace for all-day wear



CORFLEX®


Designed for life.

CORFLEXGLOBAL.COM



The CoolTex Ag Knee Wrap with Hinge provides superior comfort and support and is designed with patient-centered features that make this brace easy to apply and remove, helping to promote greater independence. Brace is available with open or closed popliteal and features bilateral heavy-duty polycentric removable hinges for added support.

The brace is manufactured from CoolTex Ag fabric to allow for easy flow of air, transfer of moisture and antimicrobial protection. Lightweight and breathable material allows for greater comfort during extended wear. Available in seven sizes from S through 4XL.



ADDITIONAL FEATURES:

- Features wrap-around design for ease of application with heavy-duty hinges for optimal medial/lateral support
- Indicated for medial/lateral instability, muscle/ligament weakness or injury, meniscal injuries and degenerative joint disease; ideal for patients with sensitivity to neoprene
- Manufactured from lightweight, breathable, antimicrobial CoolTex Ag material; open popliteal for added comfort; bilateral aluminum polycentric hinges with hyperextension stop

CoolTex™Ag 13" Anterior Closure Knee Wrap with Hinge



Open popliteal design and elastic thigh strap promote greater comfort, especially during activities that involve frequent knee flexion (also available with closed popliteal).



Brace tabs assist with application and adjustment during wear.



Colored accents make release tabs easy to see for patients with visual function difficulties.

Product Number		Size	Measurement		
Open Popliteal	Closed Popliteal		6" Above Mid-Patella	Mid-Patella	6" Below Mid-Patella
85-1303	85-1313	S	16"-18"	12"-13"	11"-13"
85-1304	85-1314	M	18"-19"	14"-16"	13"-15"
85-1305	85-1315	L	19"-21"	16"-18"	15"-17"
85-1306	85-1316	XL	21"-23"	18"-20"	17"-19"
85-1307	85-1317	2XL	23"-25"	20"-22"	19"-21"
85-1308	85-1318	3XL	25"-27"	22"-24"	21"-22"
85-1309	85-1319	4XL	27"-28"	24"-26"	22"-23"

Suggested Billing Code

L-1820 KO, elastic with condyle pads and joints



CORFLEX GLOBAL | 669 EAST INDUSTRIAL PARK DRIVE | MANCHESTER, NH 03109
CORFLEXGLOBAL.COM | TOLL FREE: 800 426 7353 | FAX: 603 623 4111