

CoolTex™ Ag 13" Knee-O-Trakker with Hinge

Silver-based fabric provides breathable comfort and antimicrobial protection



CORFLEX®

Designed for life.

CORFLEXGLOBAL.COM



The CoolTex Ag Knee-O-Trakker Brace allows the knee to move through its normal range of motion while still applying a corrective force to the patella. Dynamic C-Shaped buttress system is removable, allowing the brace to fit right or left legs, for less inventory. Textured surface helps grab the patella and is neoprene and latex-free. Brace features bilateral, aluminum hinges for added support.

The brace is manufactured from CoolTex Ag fabric to allow for easy flow of air, transfer of moisture and antimicrobial protection. Lightweight and breathable material allows for greater comfort during extended wear. Available in eight sizes from XS through 4XL.

Cool Airflow In **Heat and Perspiration Out**

Antimicrobial Silver Lining

CoolTex™ Ag
Freshness is in the Silver Lining

- A fighting layer of defense that continuously inhibits the growth of microbes, helping fabrics last longer and limiting odor formation.
- Breathable, open-mesh design promotes the transfer of heat and moisture away from the skin keeping you cool and dry.
- Nylon and spandex yarns combine to deliver superior compression and support while being extremely lightweight.

ADDITIONAL FEATURES:

- Features textured c-shaped buttress for maximum patella control and versatile strapping system that provides dynamic treatment with adjustable tension control
- Indicated for patella subluxation, patellofemoral syndrome and as part of a rehabilitation program for patella tendon injuries
- Manufactured from lightweight, breathable, antimicrobial CoolTex Ag material; bilateral, polycentric, aluminum hinges; neoprene and latex-free, EVA foam buttress

CoolTex™ Ag 13" Knee-O-Trakker with Hinge



Dynamic buttress system is contoured and textured and can be easily removed to fit left or right.



Posterior retention loop helps keep straps in position while the built-in brace tabs assist with application and adjustment.

| Product Number | Size | Measurement | | |
|----------------|------|----------------------|-------------|----------------------|
| | | 6" Above Mid-Patella | Mid-Patella | 6" Below Mid-Patella |
| 85-7272 | XS | 16"-18" | 12"-13" | 11"-12" |
| 85-7273 | S | 17"-19" | 13"-14" | 12"-13" |
| 85-7274 | M | 18"-20" | 14"-15" | 13"-14" |
| 85-7275 | L | 19"-21" | 15"-16" | 14"-15" |
| 85-7276 | XL | 20"-22" | 16"-18" | 15"-16" |
| 85-7277 | 2XL | 21"-23" | 18"-20" | 16"-17" |
| 85-7278 | 3XL | 23"-26" | 20"-22" | 17"-18" |
| 85-7279 | 4XL | 26"-29" | 22"-24" | 18"-19" |

Suggested Billing Codes

L1810 (Custom Fit)
KO, elastic with joints

L1812 (OTS)
KO, elastic with joints prefabricated



CORFLEX GLOBAL | 669 EAST INDUSTRIAL PARK DRIVE | MANCHESTER, NH 03109
CORFLEXGLOBAL.COM | TOLL FREE: 800 426 7353 | FAX: 603 623 4111