

CoolTex™ Ag 9" Wrist Splint with Strap

Antimicrobial protection combined with breathability make this the perfect splint for all-day wear



CORFLEX®

Designed for life.

CORFLEXGLOBAL.COM



The CoolTex Ag Wrist Splint is designed to provide secure and comfortable immobilization of the wrist while preserving finger mobility. Mesh liner allows splint to slip on like a glove for easy application. Palmar and dorsal stays are malleable and removable for optimal stabilization. Elastic wrist strap provides additional support and compression.

The brace is manufactured from CoolTex Ag fabric to allow for easy flow of air, transfer of moisture and antimicrobial protection. Lightweight and breathable material allows for greater comfort during extended wear. Available in seven sizes, XS through 3XL in left and right.

Cool Airflow In **Heat and Perspiration Out**

Antimicrobial Silver Lining

CoolTex™ Ag
Freshness is in the Silver Lining

- A fighting layer of defense that continuously inhibits the growth of microbes, helping fabrics last longer and limiting odor formation.
- Breathable, open-mesh design promotes the transfer of heat and moisture away from the skin keeping you cool and dry.
- Nylon and spandex yarns combine to deliver superior compression and support while being extremely lightweight.

ADDITIONAL FEATURES:

- 9" splint provides stabilization of the wrist while preserving finger mobility; elastic wrist strap provides additional compression and superior comfort
- Indicated for Carpal Tunnel Syndrome, tendonitis, wrist strains or sprains and post-cast removal
- Manufactured from lightweight, breathable, antimicrobial CoolTex Ag material; elastic; aluminum stays

Suggested Billing Code

L-3908 Wrist cock-up, non-molded

Product Number Right	Product Number Left	Size	Measurement (Wrist Circumference)
74-7600	74-7610	XS	4"-5"
74-7601	74-7611	S	5"-6"
74-7602	74-7612	M	6"-7"
74-7603	74-7613	L	7"-8"
74-7604	74-7614	XL	8"-9"
74-7605	74-7615	2XL	9"-10½"
74-7606	74-7616	3XL	10½"-12"



CORFLEX GLOBAL | 669 EAST INDUSTRIAL PARK DRIVE | MANCHESTER, NH 03109
CORFLEXGLOBAL.COM | TOLL FREE: 800 426 7353 | FAX: 603 623 4111