

# Stride OA™ Osteoarthritis Knee Brace

Prefabricated brace unloads the knee compartment helping to reduce pain associated with osteoarthritis and promote return to activity.



**CORFLEX®**

Designed for life.

CORFLEXGLOBAL.COM



Lightweight, low profile, single upright osteoarthritis knee brace features a unique angle adjustment to optimize pain relief. Laterally applied frame and hinge corrects for either medial or lateral compartment, avoids contralateral limb contact during gait and prevents tibial crest pressure issues. Perfect for bilateral brace users as well as those with limited room between the thighs or knees during gait.

Universal sizing promotes inventory reduction - only four sizes to fit all of your patients. Quick release buckles allow for easy application and removal. Compression molded liners offer comfort for extended wear. Liners and straps feature silicone grip material for optimal brace suspension.

#### ADDITIONAL FEATURES:

- Laterally applied brace helps unload the knee without interfering with gait; universal sizing for easy fitting and lower inventory requirements
- Indicated for mild to severe unicompartmental osteoarthritis (OA); unicompartmental knee conditions requiring load reduction such as articular cartilage defect repair, avascular necrosis, tibial plateau fracture or meniscus repair
- Constructed of aluminum uprights with glass-reinforced nylon cuffs; compression molded faux-suede padded liners with silicone grip

Product Number	Standard	Tall
	Patient Height up to 6'	Patient Height 6' and above
Right	81-0100	81-0102
Left	81-0101	81-0103

Fits maximum thigh circumference of 30".

Suggested Billing Code 1843 Knee Orthosis Single Upright

#### Accessories

Replacement Strap Kit	81-0090
Stride OA™ Suspension Wrap – 32" L	81-0091
Liner Kit - Right	81-0092
Liner Kit - Left	81-0093



CORFLEX GLOBAL | 669 EAST INDUSTRIAL PARK DRIVE | MANCHESTER, NH 03109  
CORFLEXGLOBAL.COM | TOLL FREE: 800 426 7353 | FAX: 603 623 4111