



**CORFLEX®**

Designed for life.

800 426 7353

## Disc Unloader Orthoses

### Lumbosacral Orthosis (LSO) & Lumbar Orthosis (LO)

*"I designed the Disc Unloader Orthosis for my patients with chronic or acute back pain. It provides relief to them by reducing pressure on the disc."*

–Vijay Vad MD

CORFLEX.COM



Disc Unloader Lumbosacral Orthosis (LSO)

Designed by a leading physiatrist, the Lumbosacral Orthosis (LSO) and Lumbar Orthosis (LO) feature a modular design that helps relieve back pain by reducing intradiscal pressure through circumferential compression. They feature a contoured rigid anterior panel for abdominal support and a rigid posterior insert set at 5° of lordosis. A posterior gel pack is also included to provide the therapeutic benefits of cold/hot therapy.

These Disc Unloader Orthoses are designed with patient comfort in mind. Made from a soft breathable lightweight fabric, the modular design of these versatile spinal orthoses promote their use throughout the continuum of care to provide stability and support without interfering with daily activities.

### Disc Unloader LSO (L0631)

Offers support from T9 vertebrae to sacrococcygeal junction

- > Orthosis comes standard with a 6" anterior panel and 15" posterior panel and cold/hot gel pack; also available with extended 10" anterior panel for additional support of pendulous abdomen
- > Constructed of soft breathable fabric with contoured rigid A.B.S. anterior and posterior plastic inserts, elastic double pulls and contact closures
- > Indicated for chronic or acute back pain, discogenic disc disease, post-surgical stabilization, spinal stenosis, herniated disc, stable non-displaced spinal fractures (L2-L5), radiculopathy, sciatica, spondylosis, spondylolisthesis, facet syndrome and general weakness of the lumbosacral spine



6" Anterior Panel



15" Posterior Insert



10" Anterior Panel

Product Number

	XS	S	M	L	XL	2XL	3XL	4XL
6" Standard	35-4200	35-4201	35-4202	35-4203	35-4204	35-4205	35-4206	35-4207
10" Extended	35-4300	35-4301	35-4302	35-4303	35-4304	35-4305	35-4306	35-4307

Sizing Guide\*

	XS	S	M	L	XL	2XL	3XL	4XL
	25"-30"	30"-35"	35"-40"	40"-45"	45"-50"	50"-55"	55"-60"	60"-65"

\*Measurement is waist circumference.

Note: Gel pack sold separately - see back

Suggested Billing Code: L-0631 Lumbar-Sacral Orthosis



ALL PRODUCTS MANUFACTURED BY CORFLEX ARE LATEX FREE.

For Care & Cleaning and Application Instructions, visit [corflex.com/discunloader](http://corflex.com/discunloader)



Disc Unloader Lumbar Orthosis (LO)

## Disc Unloader LO (L0627)

Offers support from L1 to below the L5 vertebrae

- > Orthosis comes standard with 6" anterior and 9" posterior panels; Optional 10" anterior panel is available for additional support of the pendulous abdomen; includes cold/hot gel pack
- > Constructed of soft breathable fabric with contoured rigid A.B.S. anterior and posterior plastic inserts, elastic double pulls and contact closures
- > Indicated for chronic or acute low back pain, discogenic disc disease, post surgical application, strains/sprains to lumbar muscles/ligaments and postural support

Product Number

	XS	S	M	L	XL	2XL	3XL	4XL
6" Standard	35-4000	35-4001	35-4002	35-4003	35-4004	35-4005	35-4006	35-4007
10" Extended	35-4100	35-4101	35-4102	35-4103	35-4104	35-4105	35-4106	35-4107

Sizing Guide\*

	XS	S	M	L	XL	2XL	3XL	4XL
	25"-30"	30"-35"	35"-40"	40"-45"	45"-50"	50"-55"	55"-60"	60"-65"

\*Measurement is waist circumference.  
Note: Gel pack sold separately - see below

Suggested Billing Code: L-0627 Lumbar Orthosis



6" Standard

10" Extended

9" LO  
Posterior Insert

For Care & Cleaning and Application Instructions, visit [corflex.com/discunloader](http://corflex.com/discunloader)

## Accessories



Additional Gel Pack,  
9" x 12"

Product Number

Additional Gel Pack  
9" x 12"                      72-7006



ALL PRODUCTS MANUFACTURED BY CORFLEX ARE LATEX FREE.