

# Marathon Active™

## Lace Up Ankle Stabilizer

Features the same high level of security and support you've come to trust from professional ankle taping without the costly time and expense.



**CORFLEX®**

Designed for life.



# Marathon Active™

## Lace Up Ankle Stabilizer



Long lasting, high denier nylon is ideal for sports activities or all-day wear



Customize stabilization with removable padded and contoured stays



Dual figure-8 strapping system controls inversion/eversion

- > Lace-up ankle brace features dual nylon figure-8 straps that lock the calcaneus to control abnormal inversion/eversion; ventilated fabric tongue and Achilles Comfort panel provide added comfort; wide elastic Tib/Fib strap for added proximal compression
- > Indicated for chronic ankle instability, sprains or strains and protection against inversion or eversion ankle injuries
- > Constructed of Ballistic nylon; Cooltex mesh; HDPE removable stays; polyester webbing; latex-free elastic

### Marathon Active™

#### Product Number

2XS	XS	S	M	LG	XLG	2XLG	3XLG
75-8100	75-8101	75-8102	75-8103	75-8104	75-8105	75-8106	75-8107

#### Sizing Guide Circumference

9"-10"	10"-11"	11"-12"	12"-13"	13"-14"	14"-15"	15"-16"	16"+
--------	---------	---------	---------	---------	---------	---------	------

Measurement is circumference from base of heel around. Fits either left or right foot.

#### Sizing Guide - Shoe Size

Men	6-7	7.5-9	9.5-11	11.5-13	13.5-15	15+
Narrow	2XS	XS	S	M	L	XL
Regular	XS	S	M	L	XL	2XL
Wide	S	M	L	XL	2XL	3XL
Women	5-6	6.5-8	8.5-10	10.5-12	12.5-14	14+
Narrow	2XS	XS	S	M	L	XL
Regular	XS	S	M	M	L	XL
Wide	S	M	L	L	XL	2XL

Suggested Billing Code L-1902 Ankle Foot Orthosis



# CORFLEX®

Designed for life.

800 426 7353

[CORFLEX.COM](http://CORFLEX.COM)

669 EAST INDUSTRIAL PARK DRIVE / MANCHESTER, NH 03109  
FAX 603 623 4111 / EMAIL [SERVICE@CORFLEX.COM](mailto:SERVICE@CORFLEX.COM)



ALL PRODUCTS MANUFACTURED BY CORFLEX ARE LATEX FREE.