

# Lace Align™

## Spinal TLSO

Keeping you in your comfort zone.

Our newest addition to the Lace Align family of spinal orthoses helps restrict motion in the thoracic and lumbar spine, aiding in the treatment of compression and burst fractures, kyphosis and osteoporosis. Innovative 3-point strapping system and shoulder harness provide thoracic extension and greater overall patient comfort.

Hi-tech inner lining of the TLSO incorporates Outlast®, to help manage heat and moisture for all-day comfort. Universal belt and telescoping shoulder harness minimize stocking requirements and ensure your patients can be fit in just one visit.



**CORFLEX®**

Designed for life.



Featuring Outlast® proactive heat management technology

# Lace Align™ Spinal TLSO

Keeping you in your comfort zone.



Featuring Outlast® proactive heat management technology



Dual Power-Pull Closures deliver effective compression



Belt easily adjusts to comfortably fit waists from 25"-70" in circumference.

## Lace Align™ TLSO

Product Number	UNIV
	35-3300

Sizing Guide	
	25"-70" *

*Measurement taken from waist or hip circumference (use larger measurement).  
\*Extension Kit (sold separately below) required for waist circumference over 45".*


## Accessories

Extension Kit – 16"	35-3904
Gel Pouch Only	35-3910
9" x 11" Gel Pack (each)	72-7006
9" x 11" Gel Pack (box of 10)	72-7006-010

Suggested Billing Code Thoracolumbar Sacral Orthosis, Trunk Control, Rigid Posterior Panel(s), Soft Anterior Apron L-0456 (Custom Fit), L-0457 (OTS)



Sternum strap is adjustable and removable



CoolMesh™ with



- > Posterior panel stabilizes from Scapular Spine to SC-junction; hand loops and Power-Pull Closures allow patients with limited upper extremity strength and dexterity to independently apply the brace; telescoping shoulder harness provides 13" of shoulder strap height adjustment
- > Indicated for lower back pain, sprains or strains, kyphosis secondary to osteoporosis, decompressive procedures, degenerative disc disease, bulging or herniated discs, compression fractures, post-operative stabilization, spondylolisthesis, spinal stenosis, spondylosis, or severe DJD
- > Constructed of nylon, breathable CoolMesh™ fabric, molded pull tabs, heat moldable ABS posterior panel, and hook and loop closures



# CORFLEX®

Designed for life.

800 426 7353

CORFLEX.COM

669 EAST INDUSTRIAL PARK DRIVE / MANCHESTER, NH 03109  
FAX 603 623 4111 / EMAIL SERVICE@CORFLEX.COM