

Lace Align™ Spinal Orthoses

Featuring

EnduraCool™

Transfer heat away from the skin for superior cooling.



Designed to stabilize and decompress the spine in patients recovering from surgery or injury, this innovative brace incorporates EnduraCool™, a hi-tech fabric that provides proactive temperature regulation, managing heat and moisture for all-day comfort. Removable posterior panel allows “step-down” to a traditional back brace, offering moderate support and increased patient value. Also available is an optional Gel Pouch that attaches easily to deliver cold or hot therapy following surgery, injection or injury.

The Lace Align features EnduraCool, a cutting-edge cooling technology to actively transfer heat away from the body and disperse it through the fabric. This allows for continuous cooling at key touch points, ensuring patients remain cool and comfortable.



Keeps patients
3 degrees cooler!

CORFLEX®

Designed for life.



Lace Align™ Spinal Orthoses

Indications: lower back pain, sprains or strains, spinal stenosis, spondylosis, spondylolysis, spondylolisthesis, degenerative disc disease, bulging or herniated discs, burst fractures, facet syndrome and post-operative stabilization following laminectomy, foraminotomy, IDET procedures, discectomy, fusion or decompression.

Universal Sizing

Two sizes fit waist circumference ranges of 30"-60" using the Lift-and-Slide adjustment – no trimming required

Rigid Back Panel

Removes easily for step-down modularity and can be shaped with a heat gun for added customization

Power Pull™ Closures

Durable molded pull tabs allow the patient to effortlessly apply and adjust the desired level of compression

Hand Loops

Located on both front panels to promote independent brace application and removal when grip strength is impaired

CoolMesh™ Fabric

Breathable mesh helps wick moisture and encourage airflow keeping you cool and dry

Sizing Guides

Embroidered sizing marks take the guesswork out of belt length adjustment

Integrated Anterior Panels

Shaped and slotted to help eliminate irritation when sitting or bending

Removable Side Panels

Included in the LSO Plus model for added lateral stabilization

CoolMesh™ with EnduraCool™

Proactively balances temperature swings to help prevent sweating and improve patient compliance



Spinal Orthoses

Product Number	S/M	L/XL
Lumbar Orthosis (LO)	35-3111	35-3113
Suggested Billing Code: L-0627 (Custom Fit), L-0642 (OTS)		
Lumbosacral Orthosis (LSO)	35-3201	35-3203
Suggested Billing Code: L-0631 (Custom Fit), L-0648 (OTS)		
Lumbosacral Plus Orthosis (LSO Plus)	35-3211	35-3213
Suggested Billing Code: L-0637 (Custom Fit), L-0650 (OTS)		


Sizing Guide	S/M	L/XL
	30" - 45"	40" - 60"

Measurement taken from waist or hip circumference. Use the larger of the two circumferences when selecting brace size.

Accessories

A. Replacement Back Panel (LO Model)	35-3901
B. Replacement Back Panel (LSO or LSO Plus Model)	35-3902
C. Replacement Side Panels, pair (LSO Plus Model)	35-3903
D. Lace Align Extension Kit (16")	35-3904
E. Gel Pouch with 9"x11" Gel Pack	35-3911
Gel Pouch Only	35-3910
9" x 11" Gel Pack (each)	72-7006
9" x 11" Gel Pack (box of 10)	72-7006-010

EnduraCool™



Advanced Heat Transfer

EnduraCool moves heat away from the body, providing a continuous cooling effect that adapts to your needs.

Longer-lasting Comfort

30 minutes after activity, EnduraCool proves to be more than 3° cooler than the competition.

EnduraCool elevates comfort and performance through advanced, mineral-based technologies.

